

DISCLAIMER: TO THE EXTENT ALLOWED BY LOCAL LAW, THIS INFORMATION IS PROVIDED TO YOU BY THE AMERICAN ASSOCIATION OF PHYSICISTS IN MEDICINE, A NON-PROFIT ORGANIZATION ORGANIZED TO PROMOTE THE APPLICATION OF PHYSICS TO MEDICINE AND BIOLOGY, ENCOURAGE INTEREST AND TRAINING IN MEDICAL PHYSICS AND RELATED FIELDS ("AAPM"), 'AS IS' WITHOUT WARRANTIES OR CONDITIONS OF ANY KIND, WHETHER ORAL OR WRITTEN, EXPRESS OR IMPLIED. AAPM SPECIFICALLY DISCLAIMS ANY IMPLIED WARRANTIES OR CONDITIONS OF MERCHANTABILITY, SATISFACTORY QUALITY, NONINFRINGEMENT AND FITNESS FOR A PARTICULAR PURPOSE. SOME JURISDICTIONS DO NOT ALLOW EXCLUSIONS OF IMPLIED WARRANTIES OR CONDITIONS, SO THE ABOVE EXCLUSION MAY NOT APPLY TO YOU. YOU MAY HAVE OTHER RIGHTS THAT VARY ACCORDING TO LOCAL LAW.

TO THE EXTENT ALLOWED BY LOCAL LAW, IN NO EVENT WILL AAPM OR ITS SUBSIDIARIES, AFFILIATES OR VENDORS BE LIABLE FOR DIRECT, SPECIAL, INCIDENTAL, CONSEQUENTIAL OR OTHER DAMAGES (INCLUDING LOST PROFIT, LOST DATA, OR DOWNTIME COSTS), ARISING OUT OF THE USE, INABILITY TO USE, OR THE RESULTS OF USE OF THE PROVIDED INFORMATION, WHETHER BASED IN WARRANTY, CONTRACT, TORT OR OTHER LEGAL THEORY, AND WHETHER OR NOT ADVISED OF THE POSSIBILITY OF SUCH DAMAGES. YOUR USE OF THE INFORMATION IS ENTIRELY AT YOUR OWN RISK. THIS INFORMATION IS NOT MEANT TO BE USED AS A SUBSTITUTE FOR THE REVIEW OF SCAN PROTOCOL PARAMETERS BY A QUALIFIED AND CERTIFIED PROFESSIONAL. USERS ARE CAUTIONED TO SEEK THE ADVICE OF A QUALIFIED AND CERTIFIED PROFESSIONAL BEFORE USING ANY PROTOCOL BASED ON THE PROVIDED INFORMATION. AAPM IS NOT RESPONSIBLE FOR A USER'S FAILURE TO VERIFY OR CONFIRM APPROPRIATE PERFORMANCE OF THE PROVIDED SCAN PARAMETERS. SOME JURISDICTIONS DO NOT ALLOW THE EXCLUSION OR LIMITATION OF LIABILITY FOR DAMAGES, SO THE ABOVE LIMITATION MAY NOT APPLY TO YOU.

ADULT ROUTINE CHEST CT

Indications (include but are not limited to)

- Evaluation of findings on chest radiographs or other CT exams as seen on other imaging modalities;
- Evaluation of lung and other primary thoracic malignancies, and detection and evaluation of metastatic disease;
- Evaluation for thoracic manifestations of known extrathoracic diseases;
- Evaluation of chest infections;
- Evaluation of inflammations involving the chest;
- Evaluation of abnormalities within the chest; i.e. fluid or abscess
- Evaluation of the chest wall;
- Evaluation of pleural disease;
- Evaluation of the mediastinum and lymph nodes.

Diagnostic Tasks (include but are not limited to)

- Detect nodules or masses and characterize their size and shape and relationships to organs;
- Identify abnormal aeration or expansion of the lungs;
- Detect abnormal fluid collections in the chest;
- Identify abnormal air collections both in and around the lungs;
- Detect mediastinal and paravascular masses and nodules;
- Characterize chest wall masses;
- Detect calcifications in soft tissues or the mediastinum.

Key Elements

- Contrast enhancement;
- One breath-hold (motion is problematic);
- Can reconstruct additional images for high-resolution chest CT.

Contrast

- **Oral:** None.
- **Injected:** Certain indications require administration of intravenous contrast media.
- Intravenous contrast enhancement should be performed as directed by the supervising radiologist using appropriate injection protocols and in accordance with the [ACR-SPR Practice Guideline for the Use of Intravascular Contrast Media](#) and the [ACR Manual on Contrast Media](#).

Patient Positioning

- Center the patient within the gantry; this is critical for proper functioning of AEC systems.
- Patient supine, arms above head;

Scan Range

- From top of lungs through the bottom of lungs. Instruct patient to hold breath at inspiration during entire scan.

Suspension of Respiration

- Patient should be instructed to hold his/her breath at end of inspiration.

Additional Image Reconstructions

- Certain indications may require that images be reconstructed in coronal and/or sagittal planes.
- Very thin images (approximately ≤ 1 mm) may need to be reconstructed to serve as source images for the sagittal and/or coronal reformatted images.
- Creation, use, and archival of these additional images are at the discretion of the supervising radiologist and/or departmental policy. Very large datasets may result from these additional reconstructions.

Radiation Dose Management

- AEC should be used whenever possible.
- Pay careful attention to the values selected to define the desired level of image quality (eg, Noise Index, Quality Reference mAs, Standard Deviation).
- Each manufacturer will have recommendations unique to their systems and system features. Be sure to work with your CT equipment manufacturer and a qualified medical physicist to ensure safe and appropriate operation of AEC systems.
- If more than one CT localizer radiograph is acquired, AEC systems from different manufacturers can differ with respect to which one is used to determine mA and/or kV settings. Please refer to individual manufacturer protocol instructions.

CTDI measurements and calculations

- Some manufacturers utilize a z-axis “flying focal spot”, in which two unique projections are acquired at the same z-axis table position. When this technique is used, we identify it with **. The CTDIvol on the console accurately accounts for use of this feature.

Approximate Volume CT Dose Index (CTDIvol) Values

- Approximate values for CTDIvol are listed for three different patient sizes:

| | Approx. Weight (kg) | Approx. Weight (lbs) | Approx. CTDIvol (mGy) |
|-----------------|---------------------|----------------------|-----------------------|
| Small Patient | 50-70 | 110-155 | 6-11 |
| Average Patient | 70-90 | 155-200 | 10-18 |
| Large Patient | 90-120 | 200-265 | 13-23 |

The approximate CTDIvol values are for reference only and represent a dose to the CT Dose Index phantom under very specific conditions. The CTDIvol displayed on the scanner for a patient of a given size should be similar, but not necessarily an exact match, to those listed in the above table. The provided values are all based on the 32 cm diameter “body” CTDI phantom.

It is essential that users recognize that the CTDIvol values reported on the user console prior to acquiring CT localizer radiographs on a particular patient do not represent the CTDIvol that will be delivered during that patient’s scan. CT systems rely on the CT localizer radiograph to 1) estimate the patient’s size, 2) determine the tube current settings for each tube angle and table position that will yield the requested level of image quality, and 3) calculate the average CTDIvol for the patient over the prescribed scan range. Until the CT localizer radiograph is acquired, the reported CTDIvol is not patient-specific, but is based on a generic patient size.

The CTDIvol values provided here are approximate, and are intended only to provide reference ranges for the user to consider. They are for a routine CT of an adult’s chest for the general indications given at the beginning of this document. Other indications or diagnostic tasks may have different image quality and dose requirements, and hence reasonable ranges of CTDIvol may differ according to those requirements.

In this document, a small patient is considered to be approximately 50-70 kg (110-155 lbs), an average patient approximately 70-90 kg (155-200 lbs), and a large patient 90-120 kg (200-265 lbs). However, weight is not a perfect indication of patient size. A person’s height, gender and distribution of weight across the body also must be taken into account. The thickness of the body over the area to be scanned is the best indication of patient size. Bodymass index (BMI) may also be considered:

- Underweight = BMI <18.5
- Normal weight = BMI of 18.5–24.9
- Overweight = BMI of 25–29.9
- Obesity = BMI of 30 or greater

It is recognized that the median (50th percentile) patient size for adults in the USA is larger than 70 kg. However, the 70 kg patient represents the “Reference Man”, as defined by the International Commission on Radiation Protection (ICRP), upon which AEC systems and tissue weighting factors (used for effective dose estimation) are based.

INDEX OF ADULT ROUTINE ABDOMEN-PELVIS PROTOCOLS (by manufacturer)

[GE](#)

[Hitachi](#)

[Neusoft](#)

[Philips](#)

[Siemens](#)

[Toshiba](#)

ADULT ROUTINE CHEST CT (Selected GE scanners)[\(Back to INDEX\)](#)

SCOUT: AP S60-I400; scan from top of shoulder through mid-liver, if automatic exposure control is used. PA scout if manual mA is used.

| GE (Without ASIR) | LightSpeed Ultra | BrightSpeed 16 Select | BrightSpeed 16 / Lightspeed 16 | LightSpeed Pro16 | LightSpeed VCT | Discovery CT750 HD |
|--|------------------------|-----------------------|--------------------------------|--------------------|--------------------|--------------------|
| Scan Type | Helical | Helical | Helical | Helical | Helical | Helical |
| Rotation Time (s) | 0.5 | 0.8 | 0.5 | 0.5 | 0.4 | 0.4 |
| Beam Collimation (mm) | 10 | 20 | 20 | 20 | 40 | 40 |
| Pitch | 1.35 | 1.35 | 1.35 | 1.35 | 1.375 | 1.375 |
| Speed (mm/rot) | 13.5 | 27.5 | 27.5 | 27.5 | 55.0 | 55.0 |
| kV | 120 | 120 | 120 | 120 | 120 | 120 |
| SmartmA min-max mA (ave mA) NI* | No SmartmA mA = 300 | 50 – 260 (160) | 100 – 440 (260) | 100 – 650 (270) | 100 – 650 (500) | 100 – 650 (500) |
| | | 11.57 | 11.57 | 11.57 | 13.0 | 13.0 |
| SFOV | Large | Large | Large | Large | Large | Large |

RECON 1

| | | | | | | |
|----------------|-------|-------|-------|-------|-------|-------|
| Plane | Axial | Axial | Axial | Axial | Axial | Axial |
| Algorithm | Std | Std | Std | Std | Std | Std |
| Recon Mode | Full | Full | Full | Full | Full | Full |
| Thickness (mm) | 5.0 | 5.0 | 5.0 | 5.0 | 5.0 | 5.0 |
| Interval (mm) | 5.0 | 5.0 | 5.0 | 5.0 | 5.0 | 5.0 |

RECON 2

| | | | | | | |
|----------------|-------|-------|-------|-------|-------|-------|
| Plane | Axial | Axial | Axial | Axial | Axial | Axial |
| Algorithm | Lung | Lung | Lung | Lung | Lung | Lung |
| Recon Mode | Full | Full | Full | Full | Full | Full |
| Thickness (mm) | 5.0 | 5.0 | 5.0 | 5.0 | 5.0 | 5.0 |
| Interval (mm) | 5.0 | 5.0 | 5.0 | 5.0 | 5.0 | 5.0 |

RECON 3

| | | | | | | |
|----------------|-------|-------|-------|-------|-------|-------|
| Plane | Axial | Axial | Axial | Axial | Axial | Axial |
| Algorithm | Std | Std | Std | Std | Std | Std |
| Recon Mode | Plus | Plus | Plus | Plus | Plus | Plus |
| Thickness (mm) | 1.25 | 1.25 | 1.25 | 1.25 | 1.25 | 1.25 |
| Interval (mm) | 1.25 | 1.25 | 1.25 | 1.25 | 1.25 | 1.25 |

*The Noise Index value and the primary (RECON 1) image reconstruction thickness both strongly impact CTDIvol and patient dose. See: Kanal KM et al. Impact of Operator-Selected Image Noise Index and Reconstruction Slice Thickness on Patient Radiation Dose in 64-MDCT. *AJR* 2007; 189: 219-225.

| | Approx. Weight (kg) | Approx. Weight (lbs) | Approx. CTDIvol (mGy) |
|---------------|---------------------|----------------------|-----------------------|
| Small Patient | 50-70 | 110-155 | 6-11 |
| Avg. Patient | 70-90 | 155-200 | 10-18 |
| Large Patient | 90-120 | 200-265 | 13-23 |

ADULT ROUTINE CHEST CT (selected GE scanners)(continued)**([Back to INDEX](#))**

SCOUT: AP S60-I400; scan from top of shoulder through mid-liver, if automatic exposure control is used. PA scout if manual mA is used.

| GE (With ASIR) | | VCT | Discovery CT750 HD |
|-----------------------|---------|-----------------|--------------------|
| Scan Type | | Helical | Helical |
| Rotation Time (s) | | 0.4 | 0.4 |
| Beam Collimation (mm) | | 40 | 40 |
| Pitch | | 1.375 | 1.375 |
| Speed (mm/rot) | | 55.0 | 55.0 |
| kV | | 120 | 120 |
| SmartmA | min-max | 100 – 650 (250) | 100 – 650 (250) |
| | NI* | 13.0 (DR 50%) | 16.8 |
| SFOV | | Large | Large |

RECON 1

| | | |
|----------------|-------|-------|
| Plane | Axial | Axial |
| Algorithm | Std | Std |
| Recon Mode | Full | Full |
| Thickness (mm) | 5.0 | 5.0 |
| Interval (mm) | 5.0 | 5.0 |
| ASIR | SS50 | SS50 |

RECON 2

| | | |
|----------------|-------|-------|
| Plane | Axial | Axial |
| Algorithm | Lung | Lung |
| Recon Mode | Full | Full |
| Thickness (mm) | 5.0 | 5.0 |
| Interval (mm) | 5.0 | 5.0 |
| ASIR | SS50 | SS50 |

RECON 3

| | | |
|----------------|-------|-------|
| Plane | Axial | Axial |
| Algorithm | Std | Std |
| Recon Mode | Plus | Plus |
| Thickness (mm) | 1.25 | 1.25 |
| Interval (mm) | 1.25 | 1.25 |
| ASIR | SS50 | SS50 |

*The Noise Index value and the primary (RECON 1) image reconstruction thickness both strongly impact CTDIvol and patient dose. See: Kanal KM et al. Impact of Operator-Selected Image Noise Index and Reconstruction Slice Thickness on Patient Radiation Dose in 64-MDCT. *AJR* 2007; 189: 219-225.

| | Approx. Weight (kg) | Approx. Weight (lbs) | Approx. CTDIvol (mGy) |
|---------------|---------------------|----------------------|-----------------------|
| Small Patient | 50-70 | 110-155 | 6-11 |
| Avg. Patient | 70-90 | 155-200 | 10-18 |
| Large Patient | 90-120 | 200-265 | 13-23 |

ADULT ROUTINE CHEST CT (selected HITACHI scanners)[\(Back to INDEX\)](#)**SCANOGRAM:** PA; scan from top of shoulder through mid-liver.

| HITACHI | CXR4 | ECLOS 16 | SCENARIA 64 |
|------------------------|---------|-----------|-------------|
| Scan Type | Volume | Volume | Volume |
| Rotation Time (s) | 0.8 | 0.8 | 0.75 |
| Detector Configuration | 2.5 x 4 | 1.25 x 16 | 0.625 x 64 |
| Pitch | 1.25 | 1.0625 | 0.8281 |
| Speed (mm/rot) | 6.25 | 21.25 | 33.13 |
| kV | 120 | 120 | 120 |
| mA | 225 | 225 | IntelliEC |
| Adaptive mA/IntelliEC | Yes | Yes | SD 12 |
| SFOV | 500 | 500 | 500 |
| Prep Delay (s) | 20 | 30 | 30 |

RECON 1

| Series Description | Mediastinum | Mediastinum | Mediastinum |
|--------------------|-------------|-------------|-------------|
| Type | Axial | Axial | Axial |
| Filter | 33 | 33 | 33 |
| Thickness (mm) | 5 | 5 | 5 |
| Interval (mm) | 5 | 5 | 5 |

RECON 2

| Series Description | Lung | Lung | Lung |
|--------------------|-------|-------|-------|
| Type | Axial | Axial | Axial |
| Filter | 22 | 22 | 22 |
| Thickness (mm) | 5 | 5 | 5 |
| Interval (mm) | 5 | 5 | 5 |

RECON 3

| Series Description | Thin for MPR | Thin for MPR | Thin for MPR |
|--------------------|--------------|--------------|--------------|
| Type | Axial | Axial | Axial |
| Filter | 33 | 22 | 33 |
| Thickness (mm) | 2.5 | 1.25 | 1 |
| Interval (mm) | 1.25 | 0.625 | 0.5 |

| | Approx. Weight (kg) | Approx. Weight (lbs) | Approx. CTDIvol (mGy) |
|---------------|---------------------|----------------------|-----------------------|
| Small Patient | 50-70 | 110-155 | 6-11 |
| Avg. Patient | 70-90 | 155-200 | 10-18 |
| Large Patient | 90-120 | 200-265 | 13-23 |

ADULT ABDOMEN & PELVIS ROUTINE (selected NEUSOFT scanners)[\(Back to INDEX\)](#)**SCOUT:** PA, scan from top of shoulder through mid-liver.

| NEUSOFT | Neuviz DUAL | Neuviz 16 |
|-------------------|-------------|-----------|
| Scan Type | Helical | Helical* |
| Rotation Time (s) | 1.0 | 0.6 |
| Collimation | 2 x 2.5 | 16 x .75 |
| kV | 120 | 120 |
| mA | 265 | 316 |
| Pitch | 1.50 | 0.8631 |
| FOV(mm) | 400 | 300 |

RECON 1

| | | |
|---------------------------------------|-------------------|--------|
| Series Description | Chest | Chest* |
| Image Enhancement Filter/Display Mode | Body Standard (B) | SB |
| Thickness (mm) | 5.0 | 3.0 |
| Interval (mm) | 5.0 | 3.0 |

RECON 2

| | | |
|---------------------------------------|----|--------------|
| Series Description | -- | Small RECON* |
| Image Enhancement Filter/Display Mode | -- | SB |
| Thickness (mm) | -- | 1.50 |
| Interval (mm) | -- | 0.750 |

RECON 3

| | | |
|---------------------------------------|----|-------|
| Series Description | -- | Lung* |
| Image Enhancement Filter/Display Mode | -- | EB |
| Thickness (mm) | -- | 3.00 |
| Interval (mm) | -- | 1.00 |

*ADAPTIVE FILTER is active prospectively and retrospectively.

METAL ARTIFACT REDUCTION can be activated prospectively or retrospectively (as needed).

| | Approx. Weight (kg) | Approx. Weight (lbs) | Approx. CTDIvol (mGy) |
|---------------|---------------------|----------------------|-----------------------|
| Small Patient | 50-70 | 110-155 | 6-11 |
| Avg. Patient | 70-90 | 155-200 | 10-18 |
| Large Patient | 90-120 | 200-265 | 13-23 |

ADULT ROUTINE CHEST CT (selected PHILIPS scanners)[\(Back to INDEX\)](#)**SURVIEW:** PA; scan from top of shoulder through mid-liver.

| PHILIPS | Brilliance 16 slice | Brilliance 64 slice | Ingenuity CT | Brilliance iCT SP | Brilliance iCT |
|-------------------|---------------------|---------------------|---------------|-------------------|----------------|
| Scan Type | Helical | Helical | Helical | Helical | Helical |
| Rotation Time (s) | 0.5 | 0.5 | 0.4 | 0.33 | 0.33 |
| Collimation | 16 × 1.5 mm | 64 × 0.625 mm | 64 × 0.625 mm | 64 × 0.625 mm | 128 × 0.625 mm |
| Coverage (mm) | 24 | 40 | 40 | 40 | 80 |
| kV | 120 | 120 | 120 | 120 | 120 |
| mAs | 140 | 140 | 140 | 125 | 125 |
| Pitch | 1.0 | 1.0 | 1.1 | 1.2 | 1.0 |
| FOV (mm) | 350–500 | 350–500 | 350–500 | 350–500 | 350–500 |

RECON 1

| | | | | | |
|----------------|-------|-------|-------|-------|-------|
| Type | Axial | Axial | Axial | Axial | Axial |
| Filter | B/C | B/C | B/C | B/C | B/C |
| Thickness (mm) | 3 | 3 | 3 | 3 | 3 |
| Increment (mm) | 3 | 3 | 3 | 3 | 3 |

RECON 2

| | | | | | |
|----------------|---------|---------|---------|---------|---------|
| Type | Coronal | Coronal | Coronal | Coronal | Coronal |
| Filter | B/C | B/C | B/C | B/C | B/C |
| Thickness (mm) | 3 | 3 | 3 | 3 | 3 |
| Increment (mm) | 3 | 3 | 3 | 3 | 3 |

RECON 3

| | | | | | |
|----------------|-------|-------|-------|-------|-------|
| Type | Axial | Axial | Axial | Axial | Axial |
| Filter | YA | YB | YB | YB | YB |
| Thickness (mm) | 2 | 0.9 | 0.9 | 0.9 | 0.9 |
| Increment (mm) | 1 | 0.45 | 0.45 | 0.45 | 0.45 |

RECON 4

| | | | | | |
|----------------|---------|---------|---------|---------|---------|
| Type | Coronal | Coronal | Coronal | Coronal | Coronal |
| Filter | YA | YB | YB | YB | YB |
| Thickness (mm) | 2 | 1 | 1 | 1 | 1 |
| Increment (mm) | 2 | 1 | 1 | 1 | 1 |

| | Approx. Weight (kg) | Approx. Weight (lbs) | Approx. CTDIvol (mGy) |
|---------------|---------------------|----------------------|-----------------------|
| Small Patient | 50-70 | 110-155 | 6-11 |
| Avg. Patient | 70-90 | 155-200 | 10-18 |
| Large Patient | 90-120 | 200-265 | 13-23 |

ADULT ROUTINE CHEST CT (selected SIEMENS scanners)[\(Back to INDEX\)](#)**TOPOGRAM:** PA; scan from top of shoulder through mid-liver.

| SIEMENS | Sensation 64 | Definition (Dual source 64-slice) | Definition AS+ (128-slice) | Definition Flash (Dual source 128-slice) |
|------------------------|--|--|---|---|
| Software version | VB30 | VA34 | VA44 | VA44 |
| Scan Mode | Spiral | Spiral | Spiral | Spiral |
| Rotation Time (s) | 0.5 | 0.5 | 0.5 | 0.5 |
| Detector Configuration | **64 × 0.6 mm (32 × 0.6 mm = 19.2 mm) | **64 × 0.6 mm (32 × 0.6 mm = 19.2 mm) | **128 × 0.6 mm (64 × 0.6 mm = 38.4 mm) | **128 × 0.6 mm (64 × 0.6 mm = 38.4 mm) |
| Pitch | 1.4 | 1.2 | 1.2 | 1.2 |
| kV* | 120 | 120 | 120 | 120 |
| Quality ref. mAs | 100 | 110 | 110/65*** | 110/65*** |
| CARE Dose4D | ON | ON | ON | ON |

RECON 1

| | | | | |
|----------------|-------|-------|--------------|--------------|
| Type | Axial | Axial | Axial | Axial |
| Kernel | B31f | B31f | B31f/I31f*** | B31f/I31f*** |
| Slice (mm) | 5 | 5 | 5 | 5 |
| Increment (mm) | 5 | 5 | 5 | 5 |

RECON 2

| | | | | |
|----------------|-------|-------|--------------|--------------|
| Type | Axial | Axial | Axial | Axial |
| Kernel | B70f | B70f | B70f/I70f*** | B70f/I70f*** |
| Slice (mm) | 1.5 | 1.5 | 1.5 | 1.5 |
| Increment (mm) | 1.5 | 1.5 | 1.5 | 1.5 |

RECON 3

| | | | | |
|----------------|-------|-------|--------------|--------------|
| Type | Axial | Axial | Axial | Axial |
| Kernel | B31f | B31f | B31f/I31f*** | B31f/I31f*** |
| Slice (mm) | 1.0 | 1.0 | 1.0 | 1.0 |
| Increment (mm) | 0.7 | 0.7 | 0.7 | 0.7 |

RECON 4

| | | | | |
|----------------|-------|-------|--------------|--------------|
| Type | Axial | Axial | Axial | Axial |
| Kernel | B70f | B70f | B70f/I70f*** | B70f/I70f*** |
| Slice (mm) | 1.0 | 1.0 | 1.0 | 1.0 |
| Increment (mm) | 0.7 | 0.7 | 0.7 | 0.7 |

* If scanner is equipped with automatic kV selection (CARE kV), this should be activated by selecting "On".

- For non-contrast-enhanced exams, a Reference kV of 120 and a "Dose saving optimized for" slider position of 2 is recommended.
- For contrast-enhanced exams, a Reference kV of 120 and a "Dose saving optimized for" slider position of 7 is recommended.
- For exams requiring both a non-contrast-enhanced and a contrast-enhanced scan, where the change in mean CT number of a region may be relevant to the diagnosis, it is important that the kV be the same for both scans. You can force CARE kV to use the same kV by linking the series. The optimization is then done by considering all linked acquisitions and their individual slider settings.

** Indicates that a z-axis "flying focal spot" technique is used to obtain twice as many projections per rotation as detector rows.

*** With SAFIRE.

NOTE: CARE Bolus (optional) may be used to optimize the bolus timing. Set the ROI for monitoring scan in the aorta at the level of the diaphragm with triggering threshold of 120 HU, or use manual triggering.

| | Approx. Weight (kg) | Approx. Weight (lbs) | Approx. CTDIvol (mGy) |
|---------------|---------------------|----------------------|-----------------------|
| Small Patient | 50-70 | 110-155 | 6-11 |
| Avg. Patient | 70-90 | 155-200 | 10-18 |
| Large Patient | 90-120 | 200-265 | 13-23 |

ADULT ROUTINE CHEST CT (selected TOSHIBA scanners)[\(Back to INDEX\)](#)**SCANOGRAPH:** Dual scano: PA and LAT; 500 mm above shoulders to below liver.

| TOSHIBA | Aq32 | Aq64 | AqPRIME | AqPremium | AqONE |
|------------------------|-----------------|-----------------|-----------------|-----------------|-----------------|
| Scan Type | Helical | Helical | Helical | Helical | Helical |
| Rotation Time (s) | 0.5 | 0.5 | 0.5 | 0.5 | 0.5 |
| Detector Configuration | 32 x 0.5 | 64 x 0.5 | 80 x 0.5 | 160 x 0.5 | 160 x 0.5 |
| Pitch | Fast (1.41) | Fast (1.48) | Fast (1.388) | Fast (0.994) | Fast (0.994) |
| Speed (mm/rot) | | 13.5 | | 18.75 | 39.36 |
| kV | 120 | 120 | 120 | 120 | 120 |
| mA | SURE Exposure | SURE Exposure | SURE Exposure | SURE Exposure | SURE Exposure |
| SURE Exposure | Std (SD = 12.5) | Std (SD = 12.5) | Std (SD = 12.5) | Std (SD = 12.5) | Std (SD = 12.5) |
| SFOV | LFOV | LFOV | LFOV | LFOV | LFOV |
| Breath-hold | Inspiration | Inspiration | Inspiration | Inspiration | Inspiration |
| Prep Delay (s) | 35 | 35 | 35 | 35 | 35 |

RECON 1

| | | | | | |
|-----------------|----------------|----------------|----------------|----------------|----------------|
| Type | Axial | Axial | Axial | Axial | Axial |
| SURE IQ setting | Body Std Axial | Body Std Axial | Body Std Axial | Body Std Axial | Body Std Axial |
| Thickness (mm) | 5 | 5 | 5 | 5 | 5 |
| Interval (mm) | 5 | 5 | 5 | 5 | 5 |
| DFOV (cm) | 400 | 400 | 400 | 400 | 400 |

RECON 2

| | | | | | |
|-----------------|----------------|----------------|----------------|----------------|----------------|
| Type | Axial | Axial | Axial | Axial | Axial |
| SURE IQ setting | Lung Std Axial | Lung Std Axial | Lung Std Axial | Lung Std Axial | Lung Std Axial |
| Thickness (mm) | 5 | 5 | 5 | 5 | 5 |
| Interval (mm) | 5 | 5 | 5 | 5 | 5 |
| DFOV (cm) | 400 | 400 | 400 | 400 | 400 |

VOLUME

| | | | | | |
|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|
| Type | Axial | Axial | Axial | Axial | Axial |
| SURE IQ setting | Body Std Volume | Body Std Volume | Body Std Volume | Body Std Volume | Body Std Volume |
| DFOV (cm) | 400 | 400 | 400 | 400 | 400 |

| | REFORMAT 1 | REFORMAT 2 |
|----------------|------------|------------|
| Type | Coronal | Sagittal |
| Thickness (mm) | 4 | 4 |
| Interval (mm) | 4 | 4 |

*The SURE IQ setting determines the reconstruction FC as well as other post-processing and reconstruction options, such as AIDR. The SURE IQ settings listed in this document refer to the manufacturer's default settings.

| | Approx. Weight (kg) | Approx. Weight (lbs) | Approx. CTDIvol (mGy) |
|---------------|---------------------|----------------------|-----------------------|
| Small Patient | 50-70 | 110-155 | 6-11 |
| Avg. Patient | 70-90 | 155-200 | 10-18 |
| Large Patient | 90-120 | 200-265 | 13-23 |